ITA	ENG	EE
-----	-----	----

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Prodotto

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Instruction for use

Product

Endosseous dental implant of the Integrated System for Single, Multiple, or Complete Prosthesis.

Indications

ITS dental implants can be used in completely edentulous cases in the mandible or maxilla to anchor a removable prosthesis, or fixed prosthodontics, becoming terminal intermediate abutments in bridges or single teeth.

Expected clinical benefits

The dental implants manufactured by ZIACOM ITS SRL are designed to guarantee clinical and aesthetic benefits in the field of oral implantology. In particular, the application of these medical devices allows you to prevent the loss of bone tissue, restore the patient's masticatory functions, support adjacent teeth and prevent gum disease resulting from tooth loss.

Contraindications

ITS implants cannot be used in the following cases:

- General healt conditions: cachexia, diabetes, hypertriodism, anaemia, leucopathies, hemorrhagic diathesis, osteomalacia, deforming osteitis, imperfect osteogenesis, alteration of the immune system an all systemic disease or drug therapies that compromise tissue repair such as immunosuppressors and corticosteroids. Furthermore, patients with neurological pathologies or psychotic, or with psychic instability, drug and alcohol habusers must be excluded. Heart and circulation pathologies represent a general surgical contraindications and therefore also implant therapy should be avoided. Analogously during pregnancy it is suitable to avoid surgery.
- Local conditions: Inappropriate bone quality, presence of soft tissue lesions (such as leucoplakia, lichen, stomatitis, epulis, etc.) hard tissue lesions (such as cysts, granuloma, root residues, inflammatory alterations, etc.). Inappropriate oral hygiene. Past or undergoing radiotherapy. Xerostomia. Bruxism and inappropriate occlusal conditions.
- Parent's age: dental implants must be taken into consideration only after complete bonegrowth in adolescents.

Recommendations

The surgical procedure to insert the implants is complex and highly specialized. The operator must have undergone appropriate training. Dental implants should be reserved only to patients that are sufficiently motivated and cooperative, with a good oral hygiene. Each implant site must have undergone an adequate diagnostic, clinical and radiographic evaluation. Incorrect procedures may cause implant loss and biological damages. We recommended an appropriate antibiotic coverage the surgical and post-surgical phase. The implants must be managed with tools that have been designed appositely for dental implantology and the prosthesis must be constructed using ITS components. Implant mobility, sensibility to percussion, bone loss and infections are indicators of implant failure, which must be removed.

Precautions

Complications that may follow implant surgery: ecchymosis, hemorrhage, hematoma, soft-tissue dehiscence, delayed healing, inflammation, infection, paraesthesiais, hyperesthesia, anesthesia, chronic pain due to the implant, maxillary sinus perforation, lesions to anatomical structures (nerves and vessels), alveolar atrophy in the maxilla or the mandible, oroantral or oronasal fistulas, damages to the noightboring teeth, bone fractures, implant fractures or tool fractures. Late complications may take place because of prosthetic overload, such as the fracture of the prosthetic overstructure; implant fracture, abutment loose screw and loss of osteointegration. Imperfections and peri-implantitis represent possible complications.

Attention

This device is disposable and possible re-use, in addition to being suitable for intended use, could result in sericus infections, with the possible loss of the dental implant and bone necrosis. The use of the product is exclusively limited to dentists and qualified surgeons. The ITS surgical line is subject to a continuous evolution. Not use if outer or inner packaging is open or have damage.

It is recommended to consult a doctor in case of implant mobility, pain and swelling of the surgical site.

Packaging

Sterilized by gamma rays. Disposable.

Storage

Preserve avoiding any damage to the packaging. The product fears humidity. Preserve at room temperature in a dry place. Do not re-use. The expiration date reported on the label refers to the intact product. Do not use if the package has been damaged.

Disposal

The devices are used in hospitals and must therefore be disposed of according to the laws in force regarding waste disposal and applied by the hospital. In particular, used devices which may therefore be polluted and biologically contaminated must be disposed of as special waste. Since they are intended for use in hospital facilities that normally implement this special collection, it is not necessary to provide these indications on the labelling.

Established performances

Device to anchor fixed or movable dental prosthesis.

Precautions for use

Consult an adequate panoramic X-ray and die cast models to determinate if the bone is suitable to locate the extract intervention area to individuate mandibular channels, the chin foramen, and the nose sinus cavities. Before beginning the surgical procedure, rinse the mouth with 0,2% chlorhexidine. While drilling the bone at max 1000/1500 rpm, accurately cool with a sterile saline solution. Fixture insertion must take place at 15/20 rpm, accurately cool with a torque control that is lower than 45 Ncm. Handle the fixture carefully using appropriate tools for the surgical insertion and its stability in the bone. Place the healing cap and screw carefully on the fixture's head before suturing. Prescribe to the patient a broad spectrum antibiotic cover for 5/7 days.

Attention

Place the adhesive label of the package in the patient's clinical record, to allow product traceability up to the final user.

Patient information

Follow an accurate oral hygiene at home. Take check-ups every six months or if you should have any local inconvenience where the fixture was inserted.

Side effects

There are no pharmacological side effects since there is no risk of substance release from the implant. There could be side effects such as: edema, paraesthesias, infections, and dehiscence, which are not due to the shape or material of the implant but to an improper surgical procedure.

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ENG

Information for correct functioning
Before connecting the abutment verify osteointegration carefully evaluating the following:
• absence of pain upon percussion.
• absence of movements.
• absence of radiological signs of peri-implant bone destruction.

- How to use**
1. Verify on the implant label the diameter, length and sterilization expiration date. This information must be transferred, together with the lot number, on the patient's clinical record with an adhesive label that is inside the package. Opening the external package: the non-sterile assistant provides opening the package and taking out the blister. If the blister is damaged, its content may no longer be sterile and therefore, should not be used.
 2. The non-sterile assistant tears the film letting the glass ampoule fall on the sterile cloth. From this moment, only the operator with sterile gloves will handle the device.
 3. Procedure for series CI and AI implants: the operator takes with his/her fingers the ampoule and with the other unscrews the cap. The implant is placed in the surgical site with a specific connection and pushed to the bottom of the osteotomy. Detach the transportation tool with a lateral movement. Then screw the healing cap that is located inside the package.

Instructions for Use
The device is disposable and its possible reuse, in addition to not be suitable for use, could lead to serious infections, with the loss of the implant and possible bone necrosis.

It is necessary to report any serious incident occurring in relation to the medical device supplied by us to the manufacturer and to the competent authority of the Member State in which it is based.

In accordance with art. 18 of MDR 2017/745 ZIACOM ITS SRL has provided that the patient with implants is provided with a device identification card.

“The summary relating to safety and clinical performance can be consulted on EUDAMED (<https://ec.europa.eu/Tools/eudamed>), identifying it through the basic UDI-DI codes, and until this platform is fully operational it can be requested by sending an email to the address: info@ziacom.it”.

Titanium grade 4 – ASTM F67 and ISO 5832-2:

Maximum composition limits. Percent mass fraction %	
Element	Quality 4A and 4B
Nitrogen	0,05
Carbon	0,08
Hydrogen	0,0125 *)
Iron	0,50
Oxygen	0,40
Titanium	Equilibrium
a) With the exception of billets, for which the maximum hydrogen content must be equal to 0.010% (mass fraction) and flat products for which the maximum hydrogen content must be equal to 0.015% (mass fraction).	

Titanium grade 5 Ti6Al4V – ASTM F136 and ISO 5832-3:

Maximum composition limits. Percent mass fraction %	
Element	Values
Aluminum	From 5.5 to 6.75
Vanadium	From 3,5 to 4,5
Iron	Maximum 0,3
Oxygen	Maximum 0,2
Carbon	Maximum 0,08
Nitrogen	Maximum 0,05
Hydrogen	Maximum 0,15 *)
Titanium	Equilibrium
a) With the exception of billets, for which the maximum hydrogen content must be equal to 0.010% (mass fraction).	

Symbol legend	
	Unique identification of the device
	Product code
	Lot number
	Medical device
	Manufacturing date
	Expiration date
	Sterilization by gamma rays; single sterile barrier.
	Keep dry
	Carefully read the instructions
	Do not use if the packaging is Damaged
	Single use
	Do not re-sterilize
	Attention, read the instructions

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