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Prosthetic accessory.

Indications

Teeth substitutions with cemented and screwed prostheses and anchorage of a movable prosthesis with ball attachments. The components of the implant system have been appropriately studied to be used in combination with ITS implants. For more complete details regarding the available components please refer to the ITS prosthetic catalogue.

Expected clinical benefits

The prosthetic accessories manufactured by ZIACOM ITS SRL are designed to guarantee clinical and aesthetic benefits in the field of oral implantology, in combination with the dental implants manufactured by ZIACOM ITS SRL. In particular, the application of these medical devices allows you to prevent the loss of bone tissue, restore the patient's masticatory functions, support adjacent teeth and prevent gum disease resulting from tooth loss.

Storage

Preserve avoiding damages to the package. Do not use the device if the device may pose a risk of infection.

Disposal

The devices are used in hospitals and must therefore be disposed of according to the laws in force regarding waste disposal and applied by the hospital. In particular, used devices which may therefore be polluted and biologically contaminated must be disposed of as special waste. Since they are intended for use in hospital facilities that normally implement this special collection, it is not necessary to provide these indications on the labelling.

Precautions for use

Disposable device. The reuse of the device may induce a risk of infection.

Precautions

Follow the instructions of the prosthodontic manual.

Please note that cobalt chromium is present in medical devices with code KUHG5000, KUHG5200 and UHG5000.

Attention

Before applying the abutment, accurately clean and disinfect the tool with suitable products, dry and sterilize in autoclave at 134°C at least 10 minutes. The abutments are intended to be used by specialized medical personnel. The reuse of the device may induce a risk of infection. ZIACOM ITS S.r.l. is not responsible for improper use of this device. It is recommended to consult a doctor in case of implant mobility, pain and swelling of the surgical site.

Patient information

Follow an accurate oral hygiene at home. Follow-up check-ups are recommended every six months or when any local problem occurs in the abutment implant site.

Side effects

No pharmacological side effects have been detected since the raw material that constitutes the abutments are historically inert.

Information for correct functioning

Before connecting the abutment, it is recommended to carefully verify the osteointegration by evaluating:

- absence of pain upon percussion.
- absence of movements.
- absence of radiological signs of peri-implant bone destruction.

Instructions for Use

The device is disposable and its possible reuse, in addition to not being suitable for its intended use, could lead to serious infections, with loss of the implant and possible bone necrosis.

It is necessary to report any serious incident that occurred related to the medical device supplied by us to the manufacturer and to the competent authority of the Member State in which it is located.

In accordance with art. 18 of MDR 2017/745 ZIACOM ITS SRL has provided that the patient with implants is provided with a device identification card.

"The summary relating to safety and clinical performance can be consulted on EUDAMED (<https://ec.europa.eu/Tools/eudamed>), identifying it through the basic UDI-DI codes, and until this platform is fully operational it can be requested by sending an email to the address: info@ziacom.it".

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Cobalt Chromium - CoCr28Mo ISO 5832-4:

Maximum composition limits. Percent mass fraction %	
Element	Values
Chrome	From 26,5 to 30,0
Molybdenum	From 4,5 to 7,0
Nickel	Maximum 1,0
Iron	Maximum 1,0
Carbon	Maximum 0,35
Manganese	Maximum 1,0
Hydrogen	Maximum 1,0
Cobalt	Balance

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Symbol legend

 UDI	Unique identification of the device
 REF	Product code
 LOT	Lot number
 MD	Medical Device
 M	Manufacturing date
 U	Keep dry
 I	Carefully read the instructions
 A	Attention, read the instructions
 NS	Not sterile
 X	Disposable

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